


Sterols from Olive Oil on ZB-MR2

Column: Zebtron™ ZB-MultiResidue™-2, GC Cap. Column 30 m x 0.25 mm x 0.20 µm, Ea
Phase: Proprietary Pesticides Phase
Dimensions: 30 meters x 0.25 mm x 0.2 µm
Order No: 7HG-G017-10
Oven Profile:

Carrier Gas: Constant Flow Helium, 1.3 mL/min
Injection: Splitless 0.000000000:1 0 µL @ 300°C
Detection: Refractive Index (180°C)
Analyst Note: Sample prepared using Saponification, Solid Phase Extraction (SPE), and Derivatization.

16543

**Products used in this application:****ANALYTES:**

- 1 Stigmasterol
- 2 b-Sitosterol
- 3 Betulin

