


Sterols from Olive Oil by FID on ZB-MR1**Column:** Zebtron™ ZB-MultiResidue™-1, GC Cap. Column 30 m x 0.25 mm x 0.25 µm, Ea**Phase:** Proprietary Pesticides Phase**Dimensions:** 30 meters x 0.25 mm x 0.25 µm**Order No:** 7HG-G016-11**Oven Profile:****Carrier Gas:** Constant Flow Helium, 2.5 mL/min**Injection:** Splitless 0.0000000000:1 0 µL @ 300°C**Detection:** Refractive Index (300°C)**Analyst Note:** Sample Preparation- Saponification, pH Adjustment, Solid Phase Extraction (SPE), Reconstituted, and Derivatized.**Products used in this application:**

16747

**ANALYTES:**

- 1 Campesterol
- 2 Stigmasterol
- 3 b-Sitosterol
- 4 Betulin

