


EPA 8270 Internal Standard Mix on ZB-5MSi

Column: Zebtron™ ZB-5MSi, GC Cap. Column 30 m x 0.25 mm x 0.25 µm, Ea
Phase: 5% Phenyl 95% Dimethylpolysiloxane
Dimensions: 30 meters x 0.25 mm x 0.25 µm
Order No: 7HG-G018-11
Oven Profile:

Carrier Gas: Constant Flow Helium, 0.8 mL/min
Injection: On-Column 0.0000000000:1 0 µL @ 43°C
Detection: Refractive Index (335°C)
Analyst Note: Analytes 10 ppm in Acetone

**Products used in this application:****ANALYTES:**

- 1 1,4-Dichlorobenzene-d4
- 2 Naphthalene-d8
- 3 Acenaphthene-d10
- 4 Phenanthrene-d10
- 5 Chrysene-d12
- 6 Perylene-d12

