


**FMOC-Met-OH on Lux 5 $\mu$ m Amylose-2 in RP****Column:** N/A**Dimensions:** N/A**Order No:** N/A**Elution Type:** N/A**Eluent A:** MeOH / 0.1% Formic Acid (80:20) with Formic acid (0.1%) in buffer**Gradient**      **Step No.**      **Time (min)**      **Pct A****Profile:**              **1**              0              100**Flow Rate:** 1 mL/min**Col. Temp.:** ambient**Detection:** Electrospray Mass Spec (ESMS) @ 0.0000000000 (ambient)**ANALYTES:****1** FMOC-Met-OH

Retention Time: 5.547 min

**2** FMOC-Met-OH

Retention Time: 6.959 min

