


FMOC-Ile-OH on Lux 5 μ m Amylose-2 in RP**Column:** N/A**Dimensions:** N/A**Order No:** N/A**Elution Type:** N/A**Eluent A:** ACN / 0.1% TFA (60:40) with TFA (0.1%) in buffer**Gradient****Profile:****Flow Rate:** 1 mL/min**Col. Temp.:** ambient**Detection:** Electrospray Mass Spec (ESMS) @ 0.0000000000 (ambient)

21616

**ANALYTES:**

- 1** FMOC-Ile-OH
Retention Time: 5.365 min
- 2** FMOC-Ile-OH
Retention Time: 6.7 min

