

**Amiodarone****Column:** N/A**Dimensions:** N/A**Order No:** N/A**Elution Type:** N/A**Eluent A:** 0.1% TFA in water**Eluent B:** 0.1% TFA in acetonitrile

| <b>Gradient</b> | <b>Step No.</b> | <b>Time (min)</b> | <b>Pct A</b> | <b>Pct B</b> |
|-----------------|-----------------|-------------------|--------------|--------------|
| <b>Profile:</b> | <b>1</b>        | 0                 | 95           | 5            |
|                 | <b>2</b>        | 15                | 35           | 65           |
|                 | <b>3</b>        | 30                | 5            | 95           |

**Flow Rate:** 0 mL/min**Col. Temp.:** ambient**Detection:** Electrospray Mass Spec (ESMS) @ 0.0000000000 (ambient)**Analyst Note:** SecurityGuard<sup>TM</sup> Guard Cartridge System extends column lifetime.

- SecurityGuard Cartridges, C8 4 x 3.0mm, 10/Pk Part No.: AJ0-4290

- Holder Part No.: KJ0-4282

No Chromatogram Available

**ANALYTES:**

- 1 Amiodarone I,II,III,VI,VII,VIII,X,XI,XII,XIII,XIV

