


Cold-Pressed Orange Oil on ZBWAX 60x0.25x0.25**Column:** Zebron™ ZB-WAX, GC Cap. Column 60 m x 0.25 mm x 0.25 µm, Ea**Phase:** Polyethylene Glycol (PEG)**Dimensions:** 60 meters x 0.25 mm x 0.25 µm**Order No:** 7KG-G007-11**Oven Profile:****Carrier Gas:** Constant Flow Helium, 0.8 mL/min**Injection:** Split 75.0000000000:1 0 µL @ 275°C**Detection:** Refractive Index (0°C)**Analyst Note:** Injection of Pure Orange Oil**Products used in this application:**

Cold-Pressed Orange Oil on ZBWAX 60x0.25x0.25

ANALYTES:

- 1 α -Pinene
- 2 beta-Pinene
- 3 Sabinene
- 4 3-Carene
- 5 beta-Myrcene
- 6 Limonene
- 7 beta-Phellandrene
- 8 Octanal
- 9 Limonene oxides
- 10 Citronellal
- 11 Cubebene
- 12 Decanal
- 13 Linalool
- 14 1-Octanol
- 15 Germacrene
- 16 trans-p-Mentha-2,8-dienol
- 17 cis-p-Mentha-2,8-dienol
- 18 Geranial
- 19 α -Terpineol
- 20 Dodecanal
- 21 Valencene
- 22 Citral
- 23 Carvone
- 24 Cadinene
- 25 Perillaldehyde
- 26 trans-Carveol
- 27 cis-Carveol
- 28 Perillol
- 29 Octanoic acid

