


Biodiesel Glyceride Analysis - SME on ZB-5HT**Column:** Zebtron[™] ZB-5HT, GC Cap. Column 15 m x 0.32 mm x 0.10 μ m, Ea**Phase:** 5% Phenyl 95% Dimethylpolysiloxane**Dimensions:** 15 meters x 0.32 mm x 0.1 μ m**Order No:** 7EM-G015-02**Oven Profile:****Carrier Gas:** Constant Flow Hydrogen, 11.3 mL/min**Injection:** Splitless 0.000000000:1 0 μ L @ 56°C**Detection:** Refractive Index (380°C)**Analyst Note:** Injection used PTV to 350 °C @ 1200 °C/min.
Courtesy of ASG (www.asg-analytik.de), a specialized Biofuel Laboratory!**Products used in this application:**

Biodiesel Glyceride Analysis - SME on ZB-5HT

ANALYTES:

- 1 Glycerine
- 2 Butanetriol
- 3 C16:0
- 4 C18:2
- 5 Monoolein
- 6 C18:0
- 7 Tricaprin
- 8 C34
- 9 C36
- 10 C38
- 11 C52
- 12 Triolein

