

## PEGylated IgG (N-Terminal PEG 40KDa) on BioSep3000 (2)

**Column:** BioSep™ 5 µm SEC-s3000 290 Å, LC Column 300 x 7.8 mm, Ea

**Dimensions:** 300 x 7.8 mm ID

**Order No:** 00H-2146-K0

**Elution Type:** Isocratic

**Eluent A:** 100mM Sodium Phosphate pH 6.8□□

**Gradient**      **Step No.**      **Time (min)**      **Pct A**

**Profile:**              **1**              0              100

**Flow Rate:** 1 mL/min

**Col. Temp.:** ambient


**Detection:** Laser Polarimeter @ 0.0000000000 (ambient)

**Analyst Note:** SecurityGuard™ Guard Cartridge System extends column lifetime.

- SecurityGuard Cartridges, GFC 3000 4 x 3.0mm ID, 10/Pk Part No.: AJ0-4488


- Holder Part No.: KJ0-4282

18939



### ANALYTES:

- 1 High MW PEG/IgG complex
- 2 IgG dimer / IgG / 1 PEG complex
- 3 IgG monomer
- 4 impurity
- 5 impurity
- 6 impurity



Products used in this application:

