

TCAs at neutral pH Kinetex 5u EVO C18 150x4.6mm

Column: Kinetex® 5µm EVO C18 100 Å, LC Column 150 x 4.6 mm, Ea
Dimensions: 150 x 4.6 mm ID
Order No: 00F-4633-E0
Elution Type: Isocratic
Eluent A: 20mM Sodium phosphate pH 7
Eluent B: Acetonitrile/Methanol 50/50
Gradient


Step No.	Time (min)	Pct A	Pct B
1	0	30	70

Profile:
Flow Rate: 1.25 mL/min
Col. Temp.: 30 °C
Detection: Electrospray Mass Spec (ESMS) @ 0.0000000000 (ambient)

22779



Products used in this application:



ANALYTES:

- 1 Protriptyline
- 2 Nortriptyline
- 3 Imipramine
- 4 Amitriptyline
- 5 Clomipramine

