

Ascorbic Acid in Dietary Supplement Capsules by HPLC

Column: Synergi™ 4 µm Hydro-RP 80 Å, LC Column 150 x 4.6 mm, Ea

Dimensions: 150 x 4.6 mm ID

Order No: 00F-4375-E0

Elution Type: Isocratic

Eluent A: 0.1% H₃PO₄ in Water

Gradient


Step No.	Time (min)	Pct A
----------	------------	-------

Profile:	1	0	100
-----------------	----------	----------	------------

Flow Rate: 1 mL/min

Col. Temp.: 25 °C

Detection: Electrospray Mass Spec (ESMS) @ 0.0000000000 (ambient)



ANALYTES:

1 Ascorbic Acid

Retention Time: 10 min

