

## ZB-DHA-PONA reproducibility

**Column:** Zebron ZB-DHA-PONA, GC Cap. Column 100m x 0.25mm x 0.5um , Ea

**Phase:**

**Dimensions:** 100 meters x 0.25 mm x 0.5 µm

**Order No:** 7MG-G042-17

**Oven Profile:** 150°C for 40 min

**Carrier Gas:** Constant Flow Hydrogen, 1.1 mL/min


**Injection:** Split 100:1 2.5 µL @ 250°C

**Detection:** GC/FID (325°C)

**Analyst Note:** n-Hexane



Products used in this application:



## ANALYTES:

## 1 Methane

Retention Time: 4.3 min

## 2 Undecane

Retention Time: 9.07 min

## 3 4-Chlorophenol

Retention Time: 10.76 min

## 4 Tridecane

Retention Time: 16.82 min

## 5 1-Undecanol

Retention Time: 20.66 min

## 6 Tetradecane

Retention Time: 24.5 min

©2025 Phenomenex Inc. All rights reserved.

For more information contact your Phenomenex Representative at support@phenomenex.com



## Dicyclohexylamine

Phenomenex products are available worldwide.

Retention Time: 26.33 min

## 8 Pentadecane

Retention Time: 36.81 min

[www.phenomenex.com](http://www.phenomenex.com)

support@phenomenex.com